





Product description

	 Virto M-Titanium M90/M70	 Virto M-10 NW O M90/M70/M50/M30	 Virto M-312 NW O M90/M70/M50/M30	 Virto M-312 M90/M70/M50/M30
Power level				
M	•	•	•	•
P	•	•	•	•
SP	•	•	•	•
UP**				•
Hardware configuration / options / technology				
Battery size	10	10	312	312
Push button (Program / Volume)	•	•	•	•
EasyPhone (for acoustic phone)	•	•	•	
MiniControl (for program or volume change)	•	•		
Volume control (wheel)			•	•
Telecoil	•	•	•	
Microphone technology	Omni	Omni	Omni	Directional
IP rating (without VC)	IP68 ⁸	IP68 ⁸		
Bluetooth®				•
Binaural VoiceStream Technology				•

Colors

Faceplate colors



26
Pink



22
Cocoa



28
Brown



06
Black



26
Pink



22
Cocoa



06
Black

Titanium only:
Shell color Titanium Gray & faceplate colors

Acrylic shell colors



21
Transparent



19
White



26
Pink



22
Cocoa



28
Brown



06
Black




36
Red Transparent



37
Blue Transparent

Shell styles

	 Virto M-Titanium	 Virto M-10 NW O	 Virto M-312 NW O	 Virto M-312
IIC	•	•		
CIC	•	•	•	
MC		•	•	
ITC		•	•	•
HS		•	•	•
FS		•	•	•

Performance levels

	Premium M90	Advanced M70	Standard M50	Essential M30	Features	Premium M90	Advanced M70	Standard M50	Essential M30
AutoSense OS™ 3.0					UltraZoom	•	•	•	•
Calm situation	•	•	•	•	SNR-Boost	•	•	•	
Speech in noise	•	•	•	•	DuoPhone	•	•	•	
Comfort in noise	•	•	•		SoundRecover2	•	•	•	•
Music	•	•			User Preference Tuning	•	•	•	•
Speech in loud noise*	•				Fine tuning channels	20	16	12	8
Speech in car	•				WhistleBlock	•	•	•	•
Comfort in echo	•				NoiseBlock	•	•	•	•
Media speech + mic*	•	•	•	•	WindBlock*	•	•		
Media music + mic*	•	•	•	•	EchoBlock	•			
Additional programs					SoundRelax	•	•		
Max. additional programs	3	3	3	3	QuickSync*	•	•	•	•
Comfort in echo	•				Acoustically Optimized Vent	•	•	•	•
Speech in loud noise*	•	•			Tinnitus Balance	•	•	•	•
Speech in 360°*	•	•			auto Acclimatization	•	•	•	•
Speech in noise	•	•	•	•	Environmental balance	•	•	•	•
Calm situation	•	•	•	•	RogerDirect*	•	•	•	•
Comfort in noise	•	•	•	•	Roger and directional setting*	•	•	•	
Music	•	•	•	•					
Acoustic phone	•	•	•	•					
Custom program	•	•	•	•					
Streaming programs									
RogerDirect + mic*	•	•	•	•					
PartnerMic + mic*	•	•	•	•					
Phone call + mic*	•	•	•	•					

Fitting

Software	Phonak Target 6.2 or higher
Interfaces	Noahlink Wireless, programming cable, HI-PRO, HI-PRO2

- Jansen, S. & Woodward, J. (2018). Love at first sound: the new Phonak precalculation. Phonak Insight, retrieved from www.phonakpro.com/evidence, accessed February 19th, 2020.
- Stewart, E., Rakita, L. & Drexler, J. (2019). StereoZoom Part 1: The benefit of wirelessly connected narrow directionality in Phonak hearing aids for speech intelligibility. Phonak Compendium, retrieved from <http://www.phonakpro.com/evidence>, February 19th, 2020.
- Winneke, A., Latzel, M. & Appleton-Huber (2018). Less listening- and memory effort in noisy situations with StereoZoom, Phonak Field Study News, retrieved from www.phonakpro.com/evidence, accessed February 19th, 2020.
- Rodrigues, T., & Liebe, S. (2018). AutoSense OS™ 3.0 The new & enhanced automatic operating system. Phonak Insight, retrieved from www.phonakpro.com/evidence, accessed February 19th, 2020.
- Kleine, T., Ochsner, G., Hoffman, C., Schumacher, D., Bishop, R. (2017). Biometric Calibration: improving Biometric Calibration: improving directionality for Phonak Virto™ B. Phonak Field Study News, retrieved from www.phonakpro.com/evidence, accessed February 19th, 2020.
- Bishop, R., Stewart, E., & Loyola, N. (2018). Titanium FitGuide – helping more than 50% of people get an even more discreet Virto B-Titanium. Phonak Field Study News, retrieved from www.phonakpro.com/evidence, accessed February 19th, 2020.
- Schwarzlos-Sooprayen, J.K. (2017). Deeper ear impressions with EasyView Otoblock. Phonak Field Study News, retrieved from www.phonakpro.com/evidence, accessed February 19th, 2020.
- IP68 indicates that the hearing aid is water resistant and dust tight. It survived continuous immersion in 1 meter of water for 60 minutes and 8 hours in a dust chamber as per the IEC60529 standard.

* Only for Virto M-312
 ** Only available in specific countries

The Bluetooth® word mark and logos are registered trademarks owned by the Bluetooth SIG, Inc.
 IOS is a trademark of Cisco, Inc.
 Android is a trademark of Google, Inc.